

## Operating Instructions: Infrared Viewers (IRV3 & IRX3)

### Introduction:

Power Technology's IR viewers convert otherwise invisible infrared radiation into a visible green image. These devices support both CW and pulsed sources down to picosecond–microsecond regimes without external synchronization.

**Applications (typical):** IR laser alignment; semiconductor wafer inspection (with microscope adapter); photo processing; thermal/industrial inspection; forensics; food sorting; fluid inspection; and biological illumination tasks.

### Installation:

- **Mounting (both models):** Use the ¼-20 tripod interface for hands-free operation.
- **Lens/Eyepiece:** Install the supplied C-mount objective, then set coarse focus on the objective and fine focus at the eyepiece. (Locking screws as shown in figures.)
- **Filters (NIR use):** Fit the IR cut-off/attenuation filter for near-IR viewing; when inspecting reflections, prefer a metallic diffuser. Paper tends to absorb/attenuate IR.

### Operation:

- **IRV3:** Charge via **USB-C** at **≤ 5 VDC / 500 mA**. (Do not exceed this rating.) An internal Li-ion pack powers the tube; after switching OFF, residual charge may keep the screen lit briefly—this is expected.
- **IRX3:** Install **two 3.6 V 18650 Li-ion** cells with correct polarity, close the cover securely, then power via the **ACTIVATE** button. Observe the LED battery indicator: solid = full, slow-flash = OK, fast-flash = low.


**Caution – Direct exposure:** Avoid exposing the tube face to concentrated radiation; permanent sensor damage can occur. **Important – These viewers are for indirect beam observation.** Do not use for direct beam viewing or to look for “mid-air” beams. View a diffuse target (IR card, matte metal, etc.)—over-exposure can damage the tube or reduce photocathode response.

## Operating Procedure & Control Description:

### IRX3 (Figure 1)

1. **Power:** Insert two 18650 cells, close the cover, press **ACTIVATE**.
2. **Brightness:** Use **screen brightness buttons** to set comfortable luminance without blooming.
3. **Focusing & Iris:** Rotate the **focus ring** to sharpness; adjust the **diaphragm (iris)** for contrast and depth; set **eyepiece focus** to your vision.
4. **Filters & Targets:** Use the **cut-off filter** for NIR; prefer **metallic** diffuse targets over paper.
5. **Battery Status:** **LED solid** = full; **gradual flash** = acceptable; **rapid flash** = low. Auto power-off after inactivity (see data sheet).

Figure 1: IRX3



- |   |                              |
|---|------------------------------|
| 1. LED                                  | 6. Handle                    |
| 2. Eyepiece focus ring                  | 7. Battery cover             |
| 3. Button ACTIVATE                      | 8. Tripod thread             |
| 4. Screen brightness adjustment buttons | 9. Diaphragm adjustment ring |
| 5. Filter                               | 10. Focus adjustment ring    |
|   | 11. Lens locking screw       |


## IRV3 (Figure 2)

**Power/Charge:** If needed, connect USB-C ( $\leq 5\text{ V}/500\text{ mA}$ ). Press the ON/OFF button to start. Residual screen glow after power-down is normal.

1. Focusing & Iris:
  - Rotate the objective focus ring to bring the scene to focus (1×: F/1.4 25 mm; 2×: F/1.8 50 mm).
  - Adjust the iris to balance brightness vs. depth of field; tighten the iris set screw after adjustment.
  - Fine-tune the eyepiece focus to your vision.
2. Filters & Targets: For NIR, install the IR cut-off filter. View a diffuse reflector (e.g., IR card or matte metal).
3. Tripod Use: Thread onto  $\frac{1}{4}$ -20 tripod for stable, hands-free operation.

**Notes on Visualization:** For wavelengths across the IRV model family (IRV1300/1700/2000), detection thresholds vary with wavelength and optical configuration; observing a diffuse target at 1 m typically requires the minimum irradiances indicated in the technical curves.

**Figure 2: IRV3**



- |                          |                         |
|--------------------------|-------------------------|
| 1. USB-C for charging    | 6. Handle               |
| 2. Button ON/OFF         | 7. LED                  |
| 3. Focus adjustment ring | 8. Tripod thread        |
| 4. IR Filter             | 9. Iris adjustment ring |
| 5. Eyepiece              |                         |

## **Troubleshooting (Both Models)**

- Image too dim: Open iris moderately; verify correct filter; increase target reflectivity (matte aluminum card); ensure battery charge (IRX3 LED status).
- Blooming/saturation: Stop down the iris; add neutral density attenuation for very bright sources (IRX3 accessory ND filters).
- Cannot see free-space beam: Use a card/metallic target; do not attempt to “look into” the beam path.

## **Maintenance & Service:**

These products contain no user-serviceable parts. Keep optics clean and dry. Use clean compressed air first; if needed, clean glass carefully with alcohol (isopropyl) on a lint-free swab. Avoid solvents on plastics or painted surfaces.

- Handling: Protect from mechanical shock and moisture ingress.
- Battery safety (IRX3): Follow local regulations; keep out of reach of children; never expose cells to heat/flame.

## **Warranty and Repair Return Policy:**

No return will be accepted without a factory-issued RMA number clearly marked on the package. Collect/COD returns are not accepted. For warranty returns, PTI covers outbound shipment to the customer after repair. When requesting an RMA, provide a model, serial number(s), and a concise problem description.

Opening or modifying the viewer outside documented procedures may void warranty and can require re-certification under 21 CFR 1040.10/.11.


## **Laser Safety (Read Before Use):**

- Use controls and procedures only as specified to avoid hazardous radiation exposure.
- Optical instruments increase eye hazard; never aim sources toward people; avoid eye-level beam paths.
- Always wear appropriate, wavelength-rated protective eyewear for the source under test.

## Notes on Performance & Best Practices:

- Target choice dominates visibility: Use matte aluminum/ceramic IR viewing cards or rough metal, not paper, to maximize SNR.
- Keep irradiance within the tube's comfort zone: If the scene blooms, stop down the iris and/or add ND attenuation (IRX3 accessories include ND filters).
- Free-space beams: Don't rely on scattering; use projection at 1 m and consult the model's minimum power-density vs. wavelength curves for realistic expectations.


## IRV3 Performance Charts: Spectral Sensitivity & Minimum Power Density



## IRX3 Performance Charts: Minimum Power Density & Brightness Levels

### Minimal Power density

Threshold power density dependence on wavelength. The threshold power density is defined by measuring a laser beam spot on a paper, which exhibits 20% of the overall brightness (calculated as  $255 \times 20\% = 51$ ), in contrast to the background. The measurements were taken with the viewer positioned 1.15 meters away from the piece of paper.



### Brightness levels

Normalised brightness dependence on power difference from the minimum value. The power level of 0 signifies the theoretical minimal value at which the laser beam spot becomes observable on a piece of paper. It's worth noting that the viewer exhibits lower sensitivity to laser light at 1450nm compared to 1550nm or even 1900nm.

